


Product Preparation and Recommendations



With the strength your patients need at the right price point, ZirFit® Solid Zirconia allows you to deliver beautiful fits, guaranteed. Five-axis machines mill within microns to create perfect margins, contacts and occlusions for the best possible fit and contour.


Indications

ZirFit® Solid Zirconia is indicated for posterior crowns, bridges, and screw-retained or cement-retained implant crowns. It is an esthetic alternative to PFMs or full-cast restorations. ZirFit® restorations are also ideal for patients who don't have the adequate preparation space for a PFM.


Tech tips

To ensure the correct final shade of the restoration, please make sure to define the Dentin shade on your Rx form. Dark prepped areas can have a dramatic effect on the final shade of the restoration.

Posterior Prep




Anterior Prep



Product Preparation and Recommendations

Porcelain fused to metal (PFM) crowns



Porcelain fused to metal (PFM) crowns have a long and tested track record; their ability to provide predictable, lasting results has made it a popular material selection for dentists over the years.


Indications

PFM is indicated for teeth that require both high strength and esthetics. PFM crowns are also ideal for patients who need complete coverage to block out a dark Dentin.

Tech tips

- Use a metal occlusion if patient is a bruxer.
- If you prefer a porcelain margin, use a shoulder prep.
- If the PFM is part of a precision attachment case, prep guidelines may be different. Call your lab to confirm.
- Prep guidelines indicated here are the minimum for optimal results.

Posterior Prep



Anterior Prep

